

FC2 Quick Tips

Description

- FC2 is a strong thin and flexible nitrile sheath or pouch 17 cm (6.5 inches - the same length as a male condom) in length.
- FC2 is worn by the woman during sex.
- FC2 has a flexible polyurethane ring at the closed end of the pouch (the end that is inserted into the vagina) and a soft nitrile ring at the open end of the pouch (the end that remains outside of the vagina.)
- There is silicone based lubricant on the inside of the condom, additional lubrication can be used.
- FC2 works to prevent pregnancy and STIs including HIV, by lining the vagina.

Effectiveness

- Studies show the clinical effectiveness of FC2 is equivalent or comparable to FC1.
- Over 16 years of studies have shown FC1 to be an effective barrier against unintended pregnancy and STIs, including HIV.
- FC2 is a woman-initiated method of preventing both pregnancy and STIs, including HIV.

Benefits

- FC2 can increase a woman's sense of empowerment as it is worn by the woman and offers effective dual protection.
- The man does not have to have an erection for FC2 to be used which can enhance couple intimacy.
- FC2 offers another option to the male condom and is a great choice of people with latex sensitivity.